









Hapro

secret

| | | | | | | | |
|---|---|---|---|---|--|---|---|
|  |  |  |  |  |  |  |  |
| ENG / p.2 | DEU / p.6 | NED / p.10 | FRA / p.14 | SPA / p.18 | JPN / p.22 | CHS / p.26 | C 200 |

Introduction

Congratulations on the purchase of your Hapro Seecret® anti-aging appliance! VDL Hapro B.V. produces high quality products made with care and craftsmanship. The Hapro Seecret® anti-aging appliances meet the same high standards. This new product line has been developed and assembled in accordance with strictest applicable standards and contains carefully selected, high-quality components. Based on scientifically demonstrated principles, Hapro Seecret® has been proven to vitalise and improve the healthy appearance of your skin in a relaxing manner.

Before you begin using the appliance, it is recommended that you first read this instructions and the accompanying brochure carefully. Among other things the documents contains extensive information about the Seecret® anti-aging concept, it's influence on the skin and the use of the appliance. With this information you will enjoy years of pleasure from your Hapro Seecret® anti-aging appliance.

Important

- We request you when opening the packaging to always check the product for external defects, imperfections or transport damage. Should it regrettably display defects or omissions, please contact your dealer.
- Check if the voltage indicated on the appliance corresponds to the mains voltage in your home before you connect the appliance.
- If the supply cord is damaged, it must be replaced by the manufacturer or service agent or a similarly qualified person in order to avoid a hazard.
- Water and electricity are a dangerous combination! Do not use this appliance in wet surroundings (e.g. in the bathroom or near a shower or swimming pool).
- If you have just been swimming or have just taken a shower, dry yourself thoroughly before using the appliance.
- Children should be supervised to ensure that they do not play with the appliance.
- The appliance is not to be used by children or persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction.
- Always unplug the appliance after use.
- Let the appliance cool down for approx. 15 minutes before storing it.
- Make sure the vents of the appliance remain open during use.
- Colours may fade under the influence of the sun. The same effect may occur when using this appliance.
- Protect your eyes during session using the goggles provided. Contact lenses and sun glasses are not a substitute for the goggles.
- Do not use the appliance if the timer is defective.

Using the appliance

When your appliance is properly installed and functions properly it is ready for use.

Session

As with all safe and painless skin stimulation methods, a single session will not immediately improve the skin. Aging of the skin is a process that takes place over many years, so it should be no surprise that it cannot be reversed with just one session lasting several minutes. To achieve the desired effects and see visible results a series of rejuvenation sessions is required. We recommend that you indulge in 2-3 rejuvenation sessions per week over a period of 10 to 12 weeks. Keep in mind that a 48-hour rest period is required before you begin the next session.

The results will be measurable after about one month, and you will see visible results after about three months. This differs from person to person, however, and depends greatly on the skin. After three months it may be assumed that the condition of your skin has reached a level at which the frequency of the sessions can be reduced. You can then switch to 'maintenance sessions', once every two weeks.


| Session time | Number of sessions per week | Number of weeks |
|--------------|-----------------------------|-----------------|
| 20 minutes | 2-3 | 10-12 |


- A course of treatment during a continuous 10 to 12 week period
- 2 to 3 sessions per week
- Rest period of 48 hours between sessions
- A maintenance session once every two weeks, after completing a course of treatment

As far as is known, from a scientific perspective, there is no need to limit the number of treatments with collagen light. However, based on careful study and practical tests, VDL Hapro B.V. has established these usage recommendations to help you achieve the most optimal results.

Protective goggles

The protective goggles should only be used with a Hapro Seecret appliance. They should not be used for other purposes. The protective goggles comply with protection class II according to directive 89/686/EEC.



- A** **Start**
Stand-alone mode:
 By pressing start once the “+ and –” buttons can be used to set the time. Pressing for a second time will start the unit.
Coin mode:
 The unit can be started automatically by the external timer or the coin automat. The start button can be used only when the unit has to be restarted from “pause”.
- B** **Stop**
Stand-alone mode:
 Pressing stop once will cause the unit to “pause”. The lights will go out but the ventilation will continue to run. Pressing stop twice will initiate a mandatory after-cooling phase of 3 minutes. The unit cannot be started during these 3 minutes.
Coin mode:
 Pressing stop once will cause the unit to “pause”. The lights will go out but the ventilation will continue to run. The unit will not respond to pressing stop several times. The external timer or coin automat will switch the unit off. Switching off will initiate a mandatory after-cooling phase of 3 minutes. The unit cannot be started during these 3 minutes.
-  Illustration after-cooling
- C** *Stand-alone mode:*
 a. Reduce sessiontime
 b. Reduce body cooler speed
Coin mode:
 a. Reduce body cooler speed
- D** *Stand-alone mode:*
 a. Increase sessiontime
 b. Increase body cooler speed
Coin mode:
 a. Increase body cooler speed
- E** **No function**
- F** **Fan**
 Switch the body cooler on and off

General

The unit is generally set for use without external timer (stand-alone mode) or coin automat. When the unit has to operate in combination with an external timer or a coin automat (coin mode), the internal time should be set to 0 minutes as per procedure below.

- press continuously for 10 seconds
- after 10 seconds press the “-” button and continue to press until the display indicates “00”
- press once and the unit will beep in response.

Note: - The maximum session is restricted to 30 minutes.
 - Note: Before the end of the session, the unit will give a series of beeps to alert you.

Display messages

- E1** Overrule error (Switch print and/or display print defect)
- E2** Filter glass defect (Filter glass of the face tanner is broken or glass protection is defective)
- E3** Communication error (No communication between switch print and display print)


Note: When you receive one of the above error messages on your display, you should stop using the unit and remove the plug from the socket. Contact your dealer.



This symbol will appear on the display after 500 hours. It is preferable after 500 hours to replace the lamps with new. After replacing the lamps reset the lamp hour counter.

Note: On the label on the back of the unit the lamps are stipulated that are suitable for your unit.

Resetting the hour counts



Cleaning

Always unplug the appliance and let it cool down before cleaning it.

- You can clean the outside of the appliance with a damp cloth. Make sure that no water runs into the appliance.

Never use scouring pads and abrasive liquids such as scouring agents, petrol or acetone to clean the appliance.

- Clean the protective goggles with a moist cloth.

Replacement

Lamp

The light emitted by the lamps will become less bright in time (i.e. after a few years in case of normal use). When this happens, you will notice a decrease in the capacity of the appliance. You can solve this problem by setting a slightly longer session time or by having the lamps replaced.

Always have the lamps and starters replaced by a service centre authorised by Hapro. Only they have the knowledge and skills required for this job and the original spare parts for your appliance.

Protective goggles

Additional or replacement goggles can be ordered at your local dealer.

Environment

- The lamps contain substances that may pollute the environment. When you discard the lamps, do not throw them away with the normal household waste, but dispose of them at an official collection point.
- Do not throw the appliance away with the normal household waste at the end of its life, but hand it in at an official collection point for recycling. By doing this you will help to preserve the environment.

Guarantee & service

Our appliances are manufactured with the greatest attention and are of excellent quality. Every Hapro Seecret® appliance is manufactured to meet the international IEC/EN standards. We provide a manufacturer's warranty of 1 year. Your new appliance has a manufacturer's guarantee of 1 year, which commences on the date of installation. Within the warranty period all faults that are the result of material or manufacturing errors will be repaired at no cost. Please contact your dealer for any warranty enquiries. In order to make a claim under the terms of this guarantee, you must be able to produce the original invoice or receipt bearing the date of purchase. The type and serial number of the appliance should also be indicated. The warranty is subordinate to the national regulations on product warranties.

Not covered under this warranty

Appliances which have the type, serial number and/or other identification labels removed, damaged or altered.

Lamps, starters and acrylic sheets.

Damage due to purchaser abuse, neglect, transport, improper use or improper maintenance.

Damage due to use of accessories not originally supplied by Hapro or inferior accessories.

Damage due to installation, repair or assembly not done by Hapro authorised personnel.

Damage due to fire, lightning, earthquakes or other natural disasters.

Troubleshooting

If problems should arise with this appliance and you are unable to solve them by means of the information below, please contact the nearest Hapro dealer.

When your appliance does not operate as expected this may be due to a disturbance and/or defect. Before reporting this problem to your dealer you should check whether the unit has been properly installed. The following situations apply:

- Mechanical defect. Check whether your appliance has been properly installed by once more reading the assembly instructions.
- Electrical defect. Try to correct the defect using the trouble-shooting guide below.

Note: Do not use your appliance when the ventilation is defective or is not functioning optimally or the timer is defective.

| Problem | Solution |
|-------------------------------|--|
| The appliance does not start. | <ol style="list-style-type: none"> Check that the plug is properly connected to the all socket power supply There may be a power failure. Check if the power supply works by connecting another appliance. |
| The fuse blows in use | Check the rating of the fuse. This should be: 16A slow |

| Problem | Solution |
|--|---|
| Some lamps do not burn | <ul style="list-style-type: none"> a. Check for the lamp(s) involved the starter and switch using where necessary the starter from a lamp that is functioning. b. Check the de lamp(s), by swapping those that do not work with those that do. |
| Lamps burn dimly when the unit is not in use | Reverse the switch to the power cable. |
| Poor results. | <ul style="list-style-type: none"> a. The distance between the lamps and your face is larger than prescribed b. When the lamps have been used for a long time, the capacity will decrease. This effect only becomes noticeable after several years. You can solve this problem by selecting a slightly longer session time or by having the lamps replaced. |