

# ***Harpro*** seeecret

ENG / p.2	DEU / p.7	NED / p.12	FRA / p.17	SPA / p.22	JPN / p.27	CHS / p.32	C 300
-----------	-----------	------------	------------	------------	------------	------------	-------

## Introduction

Congratulations on the purchase of your Hapro Seecret® anti-aging appliance! VDL Hapro B.V. produces high quality products made with care and craftsmanship. The Hapro Seecret® anti-aging appliances meet the same high standards. This new product line has been developed and assembled in accordance with strictest applicable standards and contains carefully selected, high-quality components. Based on scientifically demonstrated principles, Hapro Seecret® has been proven to vitalise and improve the healthy appearance of your skin in a relaxing manner.

Before you begin using the appliance, it is recommended that you first read this instructions and the accompanying brochure carefully. Among other things the documents contains extensive information about the Seecret® anti-aging concept, it's influence on the skin and the use of the appliance. With this information you will enjoy years of pleasure from your Hapro Seecret® anti-aging appliance.

## Important

The information contained in this manual regarding assembly and operation should be followed in full. Failure to do so may render the product warranty or liability on behalf of the supplier null and void. If there is any uncertainty as to whether the instructions contained in this manual have been followed properly and entirely, you should contact the supplier of the appliances.

- 1 Contact your authorised Hapro dealer immediately if the appliance does not function as outlined in this manual.
- 2 Do not operate this appliance in humid and chlorine spaces.
- 3 Always make sure the fuses can carry the load. This appliance has to be protected by ground fault circuit interrupter (GFCI) of 0.03A.
- 4 Always disconnect the power to the appliance at the wall breaker (isolator) when servicing or cleaning the in-side construction of the cabin. Maintenance for which the procedure is not described in this manual should only be carried out by a qualified electrician. Never clean the electrical components with inflammable and or chemical fluids.
- 5 Ensure that the canopy does not touch the wall when opened.
- 6 Prevent the acrylic sheets from coming into contact with sharp objects.
- 7 The appliance may not be operated if:
  - the acrylic panel in the canopy and/or bench is missing
  - one or more of the fans in the appliance are not functioning properly
  - the timing mechanism is not functioning properly
  - the temperature inside the cabin is 30°C or higher
- 8 To ensure maximum effectiveness, replace the acrylic sheet after 1200 hours.
- 9 To ensure maximum effectiveness, replace the Seecret XL 120W-R emitters after 800 hours and the Seecret 25W- R emitters after 500 hours.
- 10 Check local regulations to find out the proper procedure for lamp disposal.
- 11 DesinFact may be used to clean and disinfect the acrylic sheet. Never use cleaners that contain alcohol, as these can damage the acrylic sheet.
- 12 Emitter and other technical information of the appliance is indicated on the specification sticker.
- 13 If the supply cord is damaged, it must be replaced by the manufacturer or service agent or a similarly qualified person in order to avoid a hazard.
- 14 To replace lamps and starters, remove by turning a quarter turn to the left. Placement of new lamps and starters is done in the reverse to the above mentioned procedure.
- 15 Where applicable, clean the dust filters after every 50 hours of use. Regularly (more frequent than 50 hours) check the dust filters.
- 16 Children should be supervised to ensure that they do not play with the appliance.
- 17 The appliance is not to be used by children or persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction.
- 18 Protect your eyes during session using the goggles provided. Contact lenses and sun glasses are not a substitute for the goggles.

## Assembly instructions

Before your Hapro Seecret® appliance leaves the Hapro plants it is tested thoroughly and inspected for defects and damage. We request you when opening the packaging to always check the product for external defects, imperfections or transport damage. Should it regrettably display defects or omissions, please contact your dealer.

Your appliance should be installed carefully using these assembly instructions. You should also check whether the local power supply corresponds to that stipulated on the type panel that is stuck to the back of your unit.

Do not place your appliance in a damp area and make sure that there is adequate ventilation. When in doubt consult your dealer.

## Using the appliance

When your appliance is properly installed and functions properly it is ready for use.

### Session

As with all safe and painless skin stimulation methods, a single session will not immediately improve the skin. Aging of the skin is a process that takes place over many years, so it should be no surprise that it cannot be reversed with just one session lasting several minutes. To achieve the desired effects and see visible results a series of rejuvenation sessions is required. We recommend that you indulge in 2-3 rejuvenation sessions per week over a period of 10 to 12 weeks. Keep in mind that a 48-hour rest period is required before you begin the next session.

The results will be measurable after about one month, and you will see visible results after about three months. This differs from person to person, however, and depends greatly on the skin. After three months it may be assumed that the condition of your skin has reached a level at which the frequency of the sessions can be reduced. You can then switch to 'maintenance sessions', once every two weeks.

Session time	Number of sessions per week	Number of weeks
20 minutes	2-3	10-12


- A course of treatment during a continuous 10 to 12 week period
- 2 to 3 sessions per week
- Rest period of 48 hours between sessions
- A maintenance session once every two weeks, after completing a course of treatment

As far as is known, from a scientific perspective, there is no need to limit the number of treatments with collagen light. However, based on careful study and practical tests, VDL Hapro B.V. has established these usage recommendations to help you achieve the most optimal results.

#### Protective goggles

The protective goggles should only be used with a Hapro Secret appliance. They should not be used for other purposes. The protective goggles comply with protection class II according to directive 89/686/EEC.

#### Operation



1	Start/ Confirm	5	Accent lamps
2	Stop/ Cancel/ Pauze	6	Hidden button
3	Increase body cooling	7	Display
4	Decrease body cooling		

#### General:



Keep the "6" and the "3" buttons pressed simultaneously for 10 seconds.



After 10 seconds, the display will show "000".



Pressing the "3" button means the session time can be increased (max. 40 minutes)  
Pressing the "4" button means the session time can be decreased.



Pressing the "6" will confirm the setting.



If the controls have been set for use with an external timer or coin-operated use, the display will show an "-A-" after confirmation. If the controls have been set for use with an internal timer, the display will show an "-H-" after confirmation.



A = automatic, H = manual

*Note: If no adjustment is entered within 30 seconds, the setting is automatically confirmed. When using an external timer or a coinbox, the session time indicator must be set at 0.*

### Controls in combination with an external timer or for coin-operated use:

---



The appliance will be started automatically when the external timer or coin-operated equipment switches on. If the sunbed is set to start by using the "1" button, the appliance must be started by pressing the "1" button.



The equipment will pause when the "2" button is pressed once. If stop is pressed more than once, the equipment will remain paused. The appliance will not stop until the session time is up.


*Note: If the appliance has an external timer or is coin-operated and should be started by using the "1" button, your dealer will be able to set this for you. Control without the use of an external timer or coins:*



The session time can be decreased by pressing the "1" button once.



The maximum session time will be indicated by the flashing numbers on the display.



Use the "3" and "4" button to decrease and increase the session time.



The appliance can be started by pressing the "1" button.



The equipment will pause when the "2" button is pressed once. The lights will go out and the fan will continue for another 3 minutes. The session time will keep counting down. The session can be resumed at any time by pressing the "1" button again. When you press the "2" button for the second time, the appliance will enter the compulsory cooling stage. The session has been completed.

*Note: The session time will count down.*

### Other control functions:

---



The "3" button increases the body cooling.



The "4" button decreases the body cooling.



The "5" button turns the accent lights on and off. Error messages:

## Errors:



Error! Switch part has a defect.



Error! Communication fault between control and switch parts.



Error! Communication fault switch part.



Error! Unknown phase.

### Warning

If one of the error messages appears on the display, you must stop using the appliance and unplug it. Please contact your dealer.



To be able to read the hour counter for the lights, press the "hidden button" and "decrease body cooling" button simultaneously.



Timing indicator x 10 = working hours tubes



Press the "hidden button" to exit.



To be able to reset the hour counter for the lights, press the the "stop" and "increase body cooling" button simultaneously for 10 seconds.



Press the "hidden button" until "P90" is displayed.



Press the "increase body cooling" button until "40" is displayed.



Pressing the "start" button will confirm resetting the hour counter for the lights.


### Note






The number of hours of operation is the value on the display multiplied by 10.

## Information on Lamps

The lamps mentioned are a summary of the types that are suitable for use in the appliance.

Lamps
Seecret XL 120W-R
Seecret 25W-R



Model:	Hapro Seecret	
Type:	.....	
Serial number:	00 xx 000000	
Rated voltage:	400V/3N~/PE/50Hz	
Rated power input:	.... kW	XXXXXX
	Ballast: .....	
	Ballast: ..... W	
	Ballast: ..... W	
Made by: Hapro International B.V. PO Box 73 4420 AC Kapelle, Netherlands		

## Cleaning

Always unplug the appliance and let it cool down before cleaning it.

- You can clean the outside of the appliance with a damp cloth. Make sure that no water runs into the appliance.

Never use scouring pads and abrasive liquids such as scouring agents, petrol or acetone to clean the appliance.

- Clean the protective goggles with a moist cloth.

## Replacement

### Lamp

The light emitted by the lamps will become less bright in time (i.e. after a few years in case of normal use). When this happens, you will notice a decrease in the capacity of the appliance. You can solve this problem by setting a slightly longer session time or by having the lamps replaced.

Always have the lamps and starters replaced by a service centre authorised by Hapro. Only they have the knowledge and skills required for this job and the original spare parts for your appliance.

### Protective goggles

Additional or replacement goggles can be ordered at your local dealer.

## Environment

- The lamps contain substances that may pollute the environment. When you discard the lamps, do not throw them away with the normal household waste, but dispose of them at an official collection point.
- Do not throw the appliance away with the normal household waste at the end of its life, but hand it in at an official collection point for recycling. By doing this you will help to preserve the environment.

## Guarantee & service

Our appliances are manufactured with the greatest attention and are of excellent quality. Every Hapro Seecret® appliance is manufactured to meet the international IEC/EN standards. We provide a manufacturer's warranty of 1 year. Your new appliance has a manufacturer's guarantee of 1 year, which commences on the date of installation. Within the warranty period all faults that are the result of material or manufacturing errors will be repaired at no cost. Please contact your dealer for any warranty enquiries. In order to make a claim under the terms of this guarantee, you must be able to produce the original invoice or receipt bearing the date of purchase. The type and serial number of the appliance should also be indicated. The warranty is subordinate to the national regulations on product warranties.

### Not covered under this warranty

Appliances which have the type, serial number and/or other identification labels removed, damaged or altered.

Lamps, starters and acrylic sheets.

Damage due to purchaser abuse, neglect, transport, improper use or improper maintenance.

Damage due to use of accessories not originally supplied by Hapro or inferior accessories.

Damage due to installation, repair or assembly not done by Hapro authorised personnel.

Damage due to fire, lightning, earthquakes or other natural disasters.