














Harpro seeecret

							
ENG / p.4	DEU / p.8	NED / p.12	FRA / p.16	SPA / p.20	JPN / p.24	CHS / p.28	C 75



Introduction

Congratulations on the purchase of your Hapro Seecret® anti-aging appliance! Hapro International B.V. produces high quality products made with care and craftsmanship. The Hapro Seecret® anti-aging appliances meet the same high standards. This new product line has been developed and assembled in accordance with strictest applicable standards and contains carefully selected, high-quality components. Based on scientifically demonstrated principles, Hapro Seecret® has been proven to vitalise and improve the healthy appearance of your skin in a relaxing manner.

Before you begin using the appliance, it is recommended that you first read this instructions and the accompanying brochure carefully. Among other things the documents contains extensive information about the Seecret® anti-aging concept, it's influence on the skin and the use of the appliance. With this information you will enjoy years of pleasure from your Hapro Seecret® anti-aging appliance.

General description (Fig. 1)

- A** Storage hook for goggles
- B** Reflector
- C** Hapro Seecret Collagen lamp
- D** BCS filter
- E** Protective goggles (2 pairs)
- F** Measuring tape
- G** Mains cord with cord winding facility
- H** Fan
- I** Timer
- J** Handgrip

Important

- We request you when opening the packaging to always check the product for external defects, imperfections or transport damage. Should it regrettably display defects or omissions, please contact your dealer.
- Check if the voltage indicated on the appliance corresponds to the mains voltage in your home before you connect the appliance.
- If the supply cord is damaged, it must be replaced by the manufacturer or service agent or a similarly qualified person in order to avoid a hazard.
- Water and electricity are a dangerous combination! Do not use this appliance in wet surroundings (e.g. in the bathroom or near a shower or swimming pool).
- If you have just been swimming or have just taken a shower, dry yourself thoroughly before using the appliance.
- Children should be supervised to ensure that they do not play with the appliance.
- The appliance is not to be used by children or persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction.
- Always unplug the appliance after use.
- Let the appliance cool down for approx. 15 minutes before storing it. Wind the cord round the cord storage facility on the back of the appliance (Fig. 4).
- Make sure the vents in the back of the appliance remain open during use.

- Colours may fade under the influence of the sun. The same effect may occur when using this appliance.
- Protect your eyes during session using the goggles provided. Contact lenses and sun glasses are not a substitute for the goggles.
- Do not use the appliance if the timer is defective.
- Never use the appliance when the BCS filter is damaged, broken or missing.

Using the appliance

- 1 Place the appliance on a solid, flat and stable surface.
- 2 Fully unwind the mains cord.
- 3 Put the plug in the wall socket.
- 4 Consult the table to determine the appropriate duration of the session.
- 5 Put on the protective goggles provided.
- 6 Position yourself in front of the appliance.
- 7 Check the distance between your body and the appliance.
 - The distance between your face and the BCS filter should be 20 cm. (Fig. 2)
 - Use the measuring tape to check the distance.
- 8 Switch on the appliance by setting the timer to the required session time (Fig. 3).
 - When the set time has elapsed, the timer automatically switches off the appliance.
- 9 Always unplug the appliance after use.

Session

As with all safe and painless skin stimulation methods, a single Seecret session will not immediately improve the skin. Aging of the skin is a process that takes place over many years, so it should be no surprise that it cannot be reversed with just one session lasting several minutes. To achieve the desired effects and see visible results a series of rejuvenation sessions is required. We recommend that you indulge in 2-3 rejuvenation sessions per week over a period of 10 to 12 weeks. Keep in mind that a 48-hour rest period is required before you begin the next session.

The results will be measurable after about one month, and you will see visible results after about three months. This differs from person to person, however, and depends greatly on the skin. After three months it may be assumed that the condition of your skin has reached a level at which the frequency of the sessions can be reduced. You can then switch to 'maintenance sessions', once every two weeks.

Session time	Number of sessions per week	Number of weeks
30 minutes	2-3	10-12

- A course of treatment during a continuous 10 to 12 week period
- 2 to 3 sessions per week
- Rest period of 48 hours between sessions
- A maintenance session once every two weeks, after completing a course of treatment
-

As far as is known, from a scientific perspective, there is no need to limit the number of treatments with collagen light. However, based on careful study and practical tests,

Hapro International B.V. has established these usage recommendations to help you achieve the most optimal results.

Protective goggles

The protective goggles should only be used with a Hapro appliance. They should not be used for other purposes.

The protective goggles comply with protection class II according to directive 89/686/EEC.

Cleaning

Always unplug the appliance and let it cool down before cleaning it.

- You can clean the outside of the appliance with a damp cloth. Make sure that no water runs into the appliance.

Never use scouring pads and abrasive liquids such as scouring agents, petrol or acetone to clean the appliance.

- Clean the protective goggles with a moist cloth.
- Make the BCS filter clean with a soft cloth sprinkled with some spiritus.

Replacement

Lamp

The light emitted by the lamp will become less bright in time (i.e. after a few years in case of normal use). When this happens, you will notice a decrease in the intensity of the appliance. You can solve this problem by setting a slightly longer session time or by having the lamp replaced.

Always have the lamps and starters replaced by a service centre authorised by Hapro. Only they have the knowledge and skills required for this job and the original spare parts for your appliance.

BCS filter

Never use this appliance when the BCS filter is broken, damaged or missing. Always have the BCS filter replaced right away when it is damaged or broken.

Always have the BCS filter replaced by a service centre authorised by Hapro. Only they have the knowledge and skills required for this job and the original spare parts for your Hapro Secret appliance.

Protective goggles

Additional or replacement goggles can be ordered at your local dealer.

Environment

- The lamp of this appliance contains substances that may pollute the environment. When you discard the lamp, do not throw it away with the normal household waste, but dispose of it at an official collection point.
- Do not throw away the appliance with the normal household waste at the end of its life, but hand it in at an official collection point for recycling. By doing this, you help to preserve the environment (Fig. 5).

Guarantee & service

All defects resulting from material or production faults and which occur within the warranty period, will be repaired free of charge on production of the original sales invoice. The warranty period is 24 months. For warranty queries, please contact your dealer.

Your new Hapro Seecret C75 is intended for domestic use and not for commercial purposes. Damage as a result of commercial use, as well as tubes, starters, protective supports, are not covered by the Hapro warranty.

Malfunctioning as a result of assembly and or installation faults, improper use, lightning damage and non-original parts are excluded from the Hapro warranty.

Check carefully if the equipment has been installed properly before claiming under the warranty. In case of an unwarranted claim, you will be charged service costs.

Troubleshooting

If problems should arise with this appliance and you are unable to solve them by means of the information below, please contact the nearest Hapro dealer.

Problem	Solution
The appliance does not start.	Make sure the plug has been properly inserted into the wall socket.
	There may be a power failure. Check if the power supply works by connecting another appliance.
Poor results.	The distance between the lamps and your face is larger than prescribed.
	When the lamp has been used for a long time, the capacity will decrease. This effect only becomes noticeable after several years. You can solve this problem by selecting a slightly longer session time or by having the lamp replaced.